


Experiment no.: 3

Experiment name: Determination of urea in the blood serum.

The aim of the Experiment:

Determination of urea in the blood using Urease/Salicylate *Enzymatic colorimetric method* ENDPOINT.

Urea is hydrolyzed by urease into ammonia and carbon dioxide. The ammonia generated reacts with alkaline hypochlorite and sodium salicylate in presence of sodium nitroprusside as coupling agent to yield a green chromophore. The intensity of the color formed is proportional to the concentration of urea in the sample.



Equipment and materials used in the Experiment:

- – Photometer or colorimeter capable of measuring absorbances
- – at 600±10nm.

Constant temperature incubator set at 37 °C.

- – Cuvettes with 1-cm pathlength.
- – Pipettes to measure reagent and samples.

• REAGENT COMPOSITION

R1 Enzyme reagent. Urease > 500 U/mL. Stabilizers.

- **R2 Buffered chromogen.** Phosphate buffer 20 mmol/L pH 6.9, EDTA 2 mmol/L, sodium salicylate 60 mmol/L, sodium nitroprusside 3.4 mmol/L.

- **R3 Alkaline hypochlorite.** Sodium hypochlorite 10 mmol/L, NaOH 150 mmol/L.

- **Urea standard.** Urea 50 mg/dL (8.3 mmol/L)

Organic matrix based primary standard. Concentration value is traceable to Standard Reference Material 909b.

- **Working reagent.** Mix 1 volume of **R1** + 24 volumes of **R2**. Stable
- for 4 weeks at 2-8 °C and for 7 days at 15-25 °C.

Property of the machine:

Normal UV-Vis spectrophotometer:

Machine usage:

- Wavelength set up step.
- Blank against the solvent solution using a proper cuvette.
- Reach O.D.

Experiment procedure or protocol:

1. PROCEDURE			
1. Bring reagents and samples to room temperature.			
2. Pipette into a cuvette:			
TUBES	Blank	Sample	CAL.Standard
Working reagent	1.0 mL	1.0 mL	1.0 mL
Sample	–	10 µL	–
CAL.Standard	–	–	10 µL
3. Mix and incubate for 5 minutes at 37 °C or for 10 minutes at room temperature (16-25 °C).			
4. Pipette:			
R3	1.0 mL	1.0 mL	1.0 mL
5. Mix thoroughly and incubate the tubes for 5 minutes at 37oC or for 10 minutes at room temperature (16-25 °C).			
5. Read the absorbance (A) of the samples and the standard at 600 nm against the reagent blank.			
The color is stable for at least 2 hours protected from light.			

Experiment data and results:

Serum, plasma	
REFERENCE VALUES ⁵	
Newborns (< 10 days)	6.4 - 53.5 mg/dL (1.1 - 9.0 mmol/L)
Adults (12-60 years)	15 - 40 mg/dL (2.5 - 6.6 mmol/L)
In adults over 60 years of age, the reference interval is 17-50 mg/dL (2.8-8.3 mmol/L) and the concentrations tend to be slightly higher in males than in females.	
Urine	
Adults (normal diet)	26 - 43 g/24-h (428 - 714 mmol/24-h)
A high-protein diet causes significant increases in plasma urea concentrations and urinary excretion. It is recommended that each laboratory establishes its own reference range.	

CALCULATIONS

Serum, plasma

$A_{\text{Sample}} / A_{\text{Standard}} \times C_{\text{Standard}} = \text{mg/dL urea}$

Samples with concentrations higher than 300 mg/dL (50 mmol/L) should be diluted 1:5 with saline and assayed again. Multiply the results by 5.

Conclusion:

- Where is urea come from in the blood?
- How does it affect the health?
- What are the normal level values?
- How can you determine its quantity on the blood?
- Discuss the methodology?