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

Assay of Borax $\text{Na}_2\text{B}_4\text{O}_7 \cdot 10 \text{H}_2\text{O}$

Introduction:

Borax is sodium tetra borate, $\text{Na}_2\text{B}_4\text{O}_7 \cdot 10\text{H}_2\text{O}$. is either pure or contaminated with sodium carbonate and with free boric acid. M. wt 381.4g, Medicinally used as antiseptic.

Chemical equation and principle:

Titrate the Borax with 0.5N HCl using methyl red as indicator, boric acid is librated which is very weak acid and the indicator would not work. Mannitol is added to complex the boric acid which will increase the acidity so it can directly be titrated with NaOH. The solution should be boiled to remove CO_2 (in the precence of CO_2 the end point might be not sharp). Procedure: Weigh accurately 1 g of Borax, dissolve in 20 ml water, titrate with 0.5N HCl using 2 drops of methyl red as indicator, the color changes from yellow to pink. Record the volume required, (this volume required to convert Borax to boric acid) boil and cool the solution, add 5g of mannitol and

titrate with 1N sodium hydroxide using phenolphthalein as indicator.

Calculation:

Calculate weigh of Borax in the sample Each ml of 1N sodium hydroxide is equivalent to 0.09534g of $\text{Na}_2\text{B}_4\text{O}_7 \cdot 10\text{H}_2\text{O}$