
جامعة الانبار كلية : الصيدلة قسم : فرع الادوية والسموم اسم المادة باللغة العربية: فسلجة عملي اسم المدة باللغة الإنكليزية: . **physiology lab** المرحلة: الثانية التدريسي: مروةو شكيب ذنون عنوان المحاضرة باللغة العربية: قياس حجم الدم المضغوط وفقر الدم عنوان المحاضرة باللغة الإنكليزية:

محتوى المحاضرة

Which Cellular Elements Deficiency Is Responsible Of Anemia




Anaemia

Anaemia is the term given to the :

- 1. \downarrow number of erythrocytes or
- 2. \downarrow concentration of hemoglobin and/or
- **3.** \checkmark **Hematocrit (PCV)** as long as the total blood volume is normal.

• It can be a result of decreased erythrocytes production or rapid erythrocytes loss.



Anaemia is the term given to the :

 \downarrow number of erythrocytes or

 \downarrow concentration of hemoglobin and/or

↓ Hematocrit (PCV) as long as the total blood volume is normal.

It can be a result of decreased erythrocytes production or rapid erythrocytes loss.

Packed cell volum(PCV)

PCV is the percentage (%) of the cellular elements ; RBC, WBC & platelets in the whole blood .

Since the volume of WBC, Platelets are very less,

So the PCV value is considered equivalent to the volume of packed red cells (RBC), also called Hematocrit.


Normal value ; Male ; 42-53 %

Female 37-47%

2 Procedure according to sample collection method 1. Macro method—Wintrobe's method.

• EDTA anticoagulated whole venous blood -test tube 2. Micro method—using capillary tube (Most accurate method).

• Directly from a finger prick, to a heparin coated capillary tube



Macro method –Wintrobe s method



WINTROBE'S TUBE (FOR PCV)

Micro method – using capillary tube

Fill the **capillary tube** with blood by capillary attraction. Either from free flowing finger **punctured by a sterile lancet**/ or from a well mixed anticoagulated whole venous blood .


Seal with the **modeling clay** the empty end of the capillary tube.



Place and position the capillary tube in the **radial grooves** of **the microhematocrit centrifuge** with the sealed end away from the center (pointed toward the outside).

- Centrifuge for 5 minutes at 12000 r ,
- additional centrifugation does not pack the red blood cells more
- Express the results in percentage (%). on **Hematocrit reader**





Hematocrit reader

The hematocrit value is determined by comparing the volume of RBC's to the total volume of the whole blood sample, it is usually reported as a %. Hcrt or PCV can be used to calculate Hb concentration by ; Hb = PCV -3/ 3 Hb = PCV -2/ 3 Or Hb = PCV x 0.32 (more accurate)

- Adult males = 40% 52%.
- Adult females = 37% 47%
- Newborns= 51 60%

Children = 34 - 49%

