#### Assistant Lecturer. Zuhair A. ALrawi

#### College of Pharmacy third level


### practical lab.(5) (Lipids)

### Lipids

- It is a group of heterogeneous organic compounds that do **not dissolve** in **water** and **dissolve** with non-polar solvents such as alcohol, ether, chloroform, and Carbon tetrachloride CCl4. It is make up 50% of the composition of a living cell.
- Fats contain C, H, O, and some of them Contains P and N.

#### The most important functions of lipids

- energy source.
- Provides the body with essential fatty acids necessary for growth.
- **Protecting** the internal organs in the body such as the heart, kidneys, and spleen. keeping them from shocks.
- A source of fat-soluble vitamins (A, D, E, and K).
- It enters into the **cell structure**, especially the cell membrane.



**1- Simple Lipids :** They are esters produced from the reaction of **fatty acids** with **alcohols**.

- include:-
- Fat and Oil: It is the most common fats in nature. These are glycerol esters (triple alcohol). Hydroxyl (with monohydroxy fatty acids.
- Waxes: Such as bees wax, which are esters of longchain fatty acids with alcohols Monohydroxylate.

### **2- compound Lipids :**

- They are esters produced from the reaction of fatty acids with alcohols, and include in their composition molecular groups that contain phosphorous or nitrogen.
- Include:-
- Phospholipids.
- Glycolipids.
- Lipoproteins.

**1- Copper Acetate Test** 

This test is used to distinguish between **<u>saturated</u>** and <u>**unsaturated**</u> fatty acids and Triglycerides.

### **1- Copper Acetate Test**

### • Principle:

This experiment is based on the fact that free saturated fatty acids combine with copper acetate to form **Copper salts** In the lower water layer is **Bluish green** copper precipitate. While unsaturated fatty acids give **green** copper salts in the upper layer of petroleum ether.

### **1- Copper Acetate Test**

### Method:

- In a test tube, place a few drops of lipid and
  5 mL petroleum ether, then add 3 mL copper acetate.(CH<sub>3</sub>COO)<sub>2</sub>Cu
- Shake the mixture well, then leave the tube until the mixture separates into two layers and note the following: -


### **1- Copper Acetate Test**

#### Method:

- If the ethereal upper layer appears green and the aqueous lower layer appears blue, the fat is unsaturated fatty acid.
- If the upper layer remained clear and without color with the appearance of a greenish blue precipitate in the lower layer, the fat was saturated fatty acid.

#### **Copper Acetate Test**

### The ether layer is dissolved in the oil



the upper layer remained clear and without color

copper acetate (CH3COO)2Cus

### 2- Salkowiski reaction:

This reaction is specific to cholesterol.


#### **Principle:**

The Salkowiski's test for **cholesterol**, where it interacts with concentrated sulfuric acid (H2SO4) that absorbs water and forms two products, one dissolving in **sulfuric acid** and the other in **chloroform**.

### 2- Salkowiski's reaction:

#### Method:

- put in test tube, 2 ml of a chloroform cholesterol solution + 2 ml of slowly concentrated H2SO4 along the side of the test tube.
- Shake the tube well, then leave it until the mixture settles.
- We will notice the formation of two colored layers, the upper layer which is the red chloroform layer, and the lower the acid layer and it is yellow.
- Indicates the presence of cholesterol.



**3- Rancidity Test:** 

It is the detection of unwanted fats and oils left behind for a long time.

It is a chemical change that occurs to fats that are left over for a **long period of time**, to have a distinct taste and smell due to the fatty acids released from them.

### **3- Rancidity Test:**

### Method:-

- put in test tube, 2ml of NaOH + 3 drops of  $\alpha$ -naphthol.
- Add 2 ml of fat or any oil, then shake the mixture well.
- If rancidity is present, it gives a **red** or **pinkish** color.
- If the color appears clear, watery, there is **no rancidity**.

**Rancidity Test** No rancidity **Rancidity** 

red or pinkish color

appears clear, watery

